


G O L D S E R I E S

Radiographic Systems




*Bringing it all together
for your complete x-ray needs.*

NATIONAL HEALTHCARE DISTRIBUTORS



Gold FT system



The extended tube travel of the tubestand eliminates the need for a platform during weight-bearing knee procedures.

A wide range of cross-table procedures can be easily accomplished by rotating the tubestand on its vertical axis and utilizing the tube trunnion rings.

Economical Choices for High Performance Imaging

Gold FT Four Way Float Table System Provides Effortless Operation and Increased Productivity

- Linear bearing system "floats" the patient in four directions for quick, effortless positioning.
- Generous tabletop workspace allows for easy and comfortable patient positioning.
- The smooth, quiet lift mechanism of the optional elevating table feature provides a vertical range of 22" to 33" and permits effortless transfer of stretcher, wheelchair, casted or geriatric patients up to 400 pounds.
- Recessed foot controls within the table base offer convenient and safe hands-free operation.
- Fail-safe transverse tabletop locks provide a stable mounting surface and ensure patient safety by remaining locked even during power outages.
- Easy-access electric lock control allows quick positioning of the bucky.
- A reciprocating bucky effectively eliminates grid lines on images.

Wallstands Designed to Optimize your Imaging Capabilities

- Film receptor includes a reciprocating bucky for elimination of grid lines on images, or optional grid cabinet with fine-line grid.
- Accommodates up to a 14" x 17" cassette in either direction.
- With a precision bearing system and accurate counter-balancing, the versatile wallstand provides unencumbered vertical motion of the bucky from the floor to over six feet high. One touch of the electric lock control button secures the bucky in the desired vertical position.
- Combined with the benefits of the extended tubestand travel, the need for a platform during weight bearing knee procedures is eliminated.

Rugged Yet Versatile Tubestand Designed to Prevent Wear and Expand Procedural Capabilities

- The orthopaedic tubestand glides smoothly behind the table on precision bearings that travel on wear-resistant steel tracks. Generous track lengths are provided so that the tubestand can extend beyond the table for off-table procedures.
- Vertical travel of the tube carriage is virtually effortless as the carriage bearings coast on hard steel tracks. The range of vertical motion spans from the floor to nearly four feet over the table.
- The tube carriage also provides a transverse motion mechanism that enables the operator to position the tube ± 5 " off-center of the table midpoint. A centering detent facilitates precise repositioning of the tube.
- Easy angulation of the x-ray tube makes chest and other off-table procedures possible. The degree of rotation is clearly indicated on the angulation dial of the handle bar assembly.
- Tubestand motions – longitudinal, vertical, transverse and angulation – are controlled with electric locks that hold firmly with no drift.

Digital Imaging Ready

- NHD radiographic equipment is completely compatible with most digital imaging systems. Digital imaging technology may be added at any time. Our focus on the future secures your investment today and guarantees your investment tomorrow.



Outstanding Value in General Purpose Radiography

Gold ST Stationary Table System

The ST delivers both reliability and versatility in a compact table system. Ideal for busy clinics with demanding caseloads, the heavy-duty, all-steel design and construction easily supports patients up to 300 pounds. The optional integrated tubestand is important when space is at a premium. It's design provides continuous attachment of the tubestand to the table and installs without modification to most office walls or ceilings. A unique off-table tubestand option accommodates a full range of procedures including standing knee and ankle.

Protect Your Investment

Backed by a five-year limited parts warranty, NHD radiographic systems offer exceptional performance, reliability and long-term investment value.



NHD, Inc. is comprised of a network of regional, independently owned and operated distributors of medical imaging equipment, supplies, and services. The distributors that own NHD cover the United States with more than 100 distribution and service centers, a staff of more than 700 service technicians, and over 225 professional sales consultants and inside customer service representatives. NHD shareholder companies average more than 30 years in business in the region that they serve, maintaining an unparalleled reputation for local commitment to service. For more information visit the NHD web site at www.NHD.net.



Gold ST System



The versatile tubestand and wallstand accommodate standing knee and ankle procedures.



The Gold ST is available with an optional attached tubestand for space-saving installation without wall or ceiling modifications.



For your local NHD shareholder, please contact:
NHD Corporate Office
23250 Chagrin Blvd., Suite 430
Commerce Park 5
Beachwood, OH 44122
888.NHD.CORP
www.NHD.net