

RS 85

Tube stand without ceiling attachment



RS 85 – Tube stand without ceiling attachment

The RS 85 tube stand without ceiling attachment is easy to mount and runs on a very low floor guide. The minimum room height is 2350 mm.

All X-ray tube/collimator combinations up to a total weight of 40 kg can be accommodated to the counterbalanced tube support arm.

Central operation for all movements of the RS 85 are done at the control arm. Electromagnetic locks are provided for the longitudinal and vertical movements as well as for the rotation of the tube/collimator unit.

The easy handling of the tube stand motion speaks for itself.

The tube support arm can be rotated ± 90 degrees for lateral exposures. The unlocking is very simple and done at the control arm.

The extreme low focus-floor distance of approx. 415 mm makes exposures of the lower extremities in connection with a cassette stand or a wall bucky stand possible.

Total height of 1850 mm, tube rotation ± 180 degrees.

Longitudinal travel of the tube stand 2430 mm.

Total length of the floor guide 3012 mm.

The RS 85 tube stand without ceiling attachment is economical – versatile – reliable

